

**Zimmer Biomet® X-Ray Patient Specific Instruments for Total Knee Arthroplasty****X-PSI KNEE *Training Form***

This Training Form is intended to confirm and document surgeon understanding of the *Zimmer Biomet* Patient Specific Instruments (PSI) for the Total Knee Arthroplasty (“**X-PSI Knee**”) system, which includes *Zimmer Biomet* PSI Planner Software for X-Ray and X-PSI Mechanical Guides. By signing this form, the surgeon acknowledges having received training to prepare, and execute X-PSI case surgery using the following material (at the revision below or greater):

- X-PSI Knee System Surgical Technique (1419.1-GLBL)
- Zimmer PSI Knee Planner User Guide for X-Ray (803.307)

| | | | |
|------------------------|--|-------------|--|
| Surgeon | First name | Middle name | Last name |
| | Email address – optional | | Zimmer account number (SFDC) - optional |
| | By signing this form, the surgeon confirms that his approval or rejection of the PSI surgical plan, via electronic signature in Zimmer SMS (as detailed in page 2 of this document), is the legally binding equivalent of his handwritten signature. | | |
| | Signature | | Date (yyyy-mm-dd) |
| Sales associate | First name | Middle name | Last name |
| | Email address – optional | | Phone number – optional |
| | Distributor / territory | | I have reviewed the X-PSI Surgical Technique and Planner User guide on date: |

Referenced training material is available on SMS, ZSI or the ZimmerBiomet website. Alternatively you can contact Zimmer CAS customer support via email ZimmerPSI@ZimmerCAS.com, or phone **1-866-336-7846 (option #2)**.


Zimmer Biomet Sales Associates – please fax the completed form to ZimmerCAS Customer Service at (514) 878-3801. Alternatively, you may scan the completed form and send to zimmerpsi@zimmercas.com. In the subject line, please indicate: “*ZimmerBiomet X-PSI Surgeon Training Form*”.


Reviewing and Approving or Rejecting a PSI Surgical plan in

www.zimmersms.com

Further information can be found on the support.zimmersms.com website

1. When you log in to SMS as a surgeon, the *PSI Workflow screen* will list out cases by action required. The cases that require action are those that have surgical plans that are ready to be reviewed and approved. The corresponding PSI case status of these cases is “**Planning Ready**”.
2. Click on the case you would like to review or, if you have a specific case number, you can search for the case by entering the case number in the Case # field.
3. The link to the surgical plan is located in the Personalized Solutions section of the case. Click on the link to **Review 3D Plan** in the Planner Application or **view plan in Drive** to open the pdf file (Note that the surgical plan can only be approved through the Planner Application)



| Notes | Files | Plans | Personalized Solutions |
|---|----------------|--------------------|--------------------------------------|
| Add Product | | | |
| Left Body Side | | | |
| Product | Status | Case Details | Planner |
|  XPSI | PLANNING READY | View case in Drive | Review 3D plan View plan in Drive |

4. Once a decision has been made to approve or reject a surgical plan, click on the (**Approve / Reject**) button in the Planner application review page.
 - a. To Approve the Plan, enter any relevant comments in the comment box, select the **Approve** button. A confirmation window will appear; enter your SMS user credentials. The Planner application will be closed automatically and the case status will be updated to **Surgeon Planning Approved**.
 - b. To reject the plan, enter any relevant comments in the comment box; select the **Reject** button. A confirmation window will appear; enter your SMS user credentials. The Planner application will be closed automatically and the case status will be updated to **Surgeon Planning Rejected**.