# Zimmer® PSI Knee & Zimmer® PSI Shoulder Jigs & Bone Models Package Insert

(Valid from May 1, 2015; replaces all previous versions)

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
Before using this product, the operating surgeon should study carefully the following recommendations, warnings and instructions as well as the available product-specific information (e.g., the applicable User's Guide, product larture, written Surgical Technique). Zimmer CAS is not liable for complications arising from the use of the device outside of its indicated uses, surgical technique or judgment, product selection, and similar matters outside the control of Zimmer CAS.

1) Indications for Use: Zimmer PSI Knee & Zimmer PSI Shoulder
Zimmer PSI Knee: The Zimmer PSI Knee System is indicated as an orthopedic instrument system to assist in the positioning of knee replacement components. It involves surgical planning software used pre-operatively to plan the surgical placement of the compo nents on the basis of provided patient radiological images with identifiable placement anatomical landmarks, and surgical instrument components that include patient specific or customized guides fabricated on the basis of the surgical plan to precisely reference the placement of the implant components intra-operatively per the surgical plan.

The Zimmer PSI Knee System is to be used with the following fixed bearing knee replacement systems in accordance with their indications and contraindications: NexGen® CR, NexGen CR-Flex, NexGen CR-Flex Gender, NexGen LPS, NexGen LPS-Flex, NexGen LPS-Flex Gender, Persona® CR, and Persona PS. The patient specific guide components are intended for single-use only.

Zimmer PSI Shoulder: The Zimmer PSI Shoulder is indicated, based on patient-specific radiological images with identifiable place ment anatomical landmarks, to assist in pre-operative planning and/or intra-operative guiding of surgical instruments for shoulder replacement surgical procedures on patients not precluded from being radiologically scanned.

The Zimmer PSI Shoulder is to be used with the glenoid components of the following shoulder implant systems in accordance witl their indications and contraindications: Zimmer® Bigliani/Flatow® Complete Shoulder Solution, Zimmer® Trabecular Metal™ Glenoid, and Zimmer® Trabecular Metal™ Reverse Shoulder. The Zimmer PSI Shoulder instrument guides and bone model are intended for single use only.

he Zimmer PSI Shoulder should not be used when the patient has metallic devices implanted that could interfere with the CT scan quality. Additionally, the Zimmer PSI Shoulder should not be used in cases where native bone is absent, or where a custom bone augment/graft will be used, on surfaces intended to mate with the Zimmer PSI Shoulder instrument guide:

The following are general precautions and warnings related to the use of Zimmer PSI Knee and Zimmer PSI Shoulder Jigs and Bone Models from Zimmer CAS, referred to commonly as PSI Components within this insert. Other specific warnings and precautions are included in the Surgical Technique or User Guides.

- · PSI Components should only be used by trained and qualified surgeons. The surgeon should follow a training course given by
- Zimmer or the distributor prior to use.
  PSI Components are to be used with the given implant system per the related pre-operative planning. The implant must be used in accordance with its respective package labeling. The user should refer to the surgical technique published by the implant manu-
- PSI Components are provided non-sterile. The components must be cleaned and sterilized before the surgery.
   PSI Components are single use only and they should be discarded after surgery.
- · PSI Components cannot be re-sterilized and re-cleaned if they have been soiled in any way including contact with blood, protein or
- any contaminants. They can only be re-sterilized and re-cleaned once per the cleaning and sterilization methods below only if they have not been soiled as described. This is to avoid patient infection or contamination.

  • PSI Components are designed and are specific to the patient's anatomy as it was at the moment when the patient radiological im-
- ages were acquired. If the anatomy or condition of the articular surfaces changed since the radiological images were acquired, the PSI Components should not be used.
- PSI Components instruments should not be used after the expiration date as indicated on the package label PSI Components should not be used if they have been dropped on the floor.
- . The user must verify that the PSI Components are in good condition prior to the surgery, including that they are free of bends, dents,
- PSI Components should be handled carefully and they should not be placed under other heavy instruments or equipment to avoid
- Minimize excessive heat buildup due to friction between metallic components and adjacent surfaces, whether it is bone or other
- metallic surfaces. Excess heat buildup in instruments in contact with bone may cause heat necrosis and cause deformation of the
- PSI Components should not be used if the unique patient identifier (PSI Case ID) is not legible on the Jigs.
- PSI Components that have been used and are being returned to Zimmer must be cleaned and sterilized prior to return, and be accompanied by a written confirmation of sterilization.

- 4) CleaningPSI Components should not be cleaned with metal brushes or other cleaning equipment that can damage its surface. Use a soft
- Use a water jet to flush difficult to access areas and closely mated surfaces.
- PSI Components must be sonicated (ultrasound cleaning) in an enzyme solution with a minimum cycle times of 5 minutes.
   PSI Components must be sonicated (ultrasound cleaning) in an enzyme solution with a minimum cycle times of 5 minutes.
   PSI Components must not come into contact with substances containing chlorine or fluorine. Instruments that are made wholly or partly of plastic must not come into contact with strong acids or lyes, organic or ammonia-containing solvents, aromatic and/or halogen hydrocarbons or oxidizing chemicals.

PSI Components require the Steam Sterilization method as follows:

PSI Components should not be sterilized in the packaging supplied with them. They must be cleaned prior to sterilization. The packaging method must be suitable for the sterilization process. All sterilizations should be performed using standard and regularly

### Steam Sterilization (Autoclave)

Cycle Type	Temperature <sup>1</sup>	Exposure Time <sup>1</sup>	Minimum Dry Time <sup>2</sup>	Minimum Cool Time <sup>3</sup>
Pre-Vacuum	132 °C (270 °F)	4 minutes	30 minutes	30 minutes

<sup>1</sup> Both the given cycle temperature and time can be increased to 134°C + 3°C (273.2°F + 5.4°F) and 18 minutes according to local requirements outside of the United States such as in the European Union.

<sup>2</sup> Drying times vary according to load size and should be increased for larger loads.

Cooling times vary according to the type of sterilizer used, device design, temperature and humidity of ambient environment, and type of packaging used. Cooling process should comply with ANSI/AAMI ST79.

0413 The CE mark is valid only if it is printed on the package label

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