


Nuance[®] PowerShare[™] Network

Expedite your Transfer with PowerShare

Utilize your account access to immediately share images and reports.




1

Go to www1.nuancepowershare.com and Login.

NOTE: If you do not know your password, select **Forgot Your Password?** for reset instructions.

NOTE: If this is your first time logging in, you will need to click **Accept** to accept the Terms of Use.




2

After logging in, click **Upload Images**.

When you are prompted to run the **PowerShare Java Web Start Uploader** application, click **Run**.

NOTE: You must have Java version 1.7.21 or higher installed on your computer to upload DICOM images. Java 8 Update 91 is recommended for optimal performance.



3

A new Uploader window will display. From the dropdown list, select the account you would like to upload to : **Receiving Site Name**
And upload mode : **Batch Upload**


...continued next page

Questions? Call 1-866-809-4746

Nuance® PowerShare™ Network

Expedite your Transfer with PowerShare

Utilize your account access to immediately share images and reports.




4

Select the specific study or studies you wish to upload by clicking on them. You can upload DICOM, Non-DICOM or a combination of both image types in a single upload. Non-DICOM images can be found on the **Other Images Tab**.

NOTE: Red indicates that the study is missing required information and must be edited.

NOTE: Non-DICOM images will require you to select the image and then the **Add Image Details** button to add required patient information.



5

Your upload progress will be displayed. You may now select to **Upload Another** or **Done**.

Once complete, both you and the requesting site will be able to view the images.



We will work with your site to configure a PowerShare Accelerator

enabling you to send directly from your PACS and save valuable time on your next transfer.

Questions? Call 1-866-809-4746